

ULTRA - GEL AQUARIUS 101

is used by governmental hospitals, institutions, public & private health services.

Directions: Apply to contact area, transducer or electrodes,
Remove with a tissue upon completion

Ingredients: Carbomer, Triethanolamine, Preservative, Water.

For Ultrasound Procedures Physiotherapy Electrodes

High-Transmission, Multi-Purpose
Gentle, Water-soluble, Colorless
Salt & Oil Free,
Non Irritating
Used Worldwide

Manufactured by Dakar Medical Ltd.
under a quality system certified as complying
with Iso 9001: 2008, Iso 13485: 2003
by the standards institution of Israel.
Ministry of Health License



350 ml. (11.8 oz)

External use only

